

1 Objectives

Create agar phantom for laser heating and MRTI.

2 Work

2.1 Original Phantom Recipe

Egg White Phantom:

- 1.5% agar
 - 15 g agar for 1-L phantom volume (can go as low as 10 g, but not as firm), i.e. 1.5 g per 100 mL
 - food grade agar powder ~\$8 from Amazon
 - clinical grade (white bottle, red cap) ~\$250
- 2/3 volume distilled water
- 1/3 volume liquid egg whites

Directions:

- Keep egg whites warm (under 60°C) - hot plate with stir bar in beaker ~ 50°C
- Put agar powder in water before heating - hot plate with stir bar in beaker
- Heat water and agar to 80-90°C to dissolve agar powder, will become translucent when ready (can use 2510 Branson sonicare to eliminate bubbles while stirring)
- Let water and agar cool to ~ 60°C while stirring
- Pour egg whites into beaker with water and agar
- Remove stir bars from beaker
- Slowly pour solution into container and let cool until it sets

2.2 Second Scan Recipes

1 Agar Phantom:

- 1.5% agar ($T_1 = 2640$ ms) (Mathur-De Vre et al, 1985)
- 350 mL total volume \rightarrow 5.25 g agar
- cylindrical hole left for fat reference (allow to set)
- trimmed container top first, covered with Al foil

1 Agar + Gd Phantom:

- 1.5% agar and 0.14 mM Gd
 - $T_1 = 791$ ms (Walker et al, 1989)
 - (White matter $T_1 = 832$ ms) (Wansapura et al, 1999)
- 350 mL total volume \rightarrow 5.25 g agar
- Gd dilution:
 - $C_1 V_1 = C_2 V_2$
 - $C_1 = 2.43$ mM Gd

$$- V_1 = C_2 V_2 / C_1 = (0.14 \text{ mM Gd}) (350 \text{ mL}) / (2.43 \text{ mM Gd}) = 20.16 \text{ mL Gadodiamide}$$

- cylindrical hole left for fat reference (allow to set)
- trimmed container top first, covered with Al foil

Fat References:

- 1.5% agar
- 200 mL total volume \longrightarrow 3.0 g agar
- 1:1 water:Crisco \longrightarrow 100 mL water, 100 mL Crisco
- 1 Lecithin pill to dissolve Crisco in water
- heat and stir until 80°C and complete dissolution
- pour slowly and fill cylindrical hole left in phantom

3 Directions

Heat phantom and scan to acquire full k-space data in time.